

Aurum[®] Plus

Product range



The
POSITIVE
choice

A range of stoma pouches
designed around what you value
– choice, independence, healthy
skin & leak free days.

Aurum[®] Plus 



The
TRUSTED
choice

Designed to reduce the risk
of leak channels forming.



The
STYLISH
choice

Choose between
pouches in black or sand.

There's more than one plus...



The
SKIN-FRIENDLY
choice


A hydrocolloid flange
with Manuka honey.



The
CLINICAL
choice

Available in flat and soft convex pouches,
with a choice of closed, drainable and urostomy.

Designed for people who know that the difference is in the detail.

 Water repellent pouch material

 Hook-to-hook fastening

The POSITIVE choice


Aurum[®] Plus 



Dual-Carb[®] i filter (Drainable pouch)
Dual-Carb[®] 2 filter (Closed pouch)

 6mm soft convex flange  Soft belt loops


Hydrocolloid flange with Manuka honey 



 Designed to reduce the risk of leak channels forming



 Split back cover

Tap pocket 

 New pouch profile reduces the risk of the outlet twisting and blocking the flow of urine





Aurum® Plus Soft Convex Closed

KEY

-  Black colour pouch
-  Sand colour pouch
-  Clear pouch



Product Code	Pack Size	Pouch Size	Plateau Size	Aperture Size
X1C1N126113	10	Midi	25mm	13mm Starter Hole
X1C1N126120	10	Midi	25mm	20mm
X1C1N126125	10	Midi	25mm	25mm
X1C1N226113	10	Midi	25mm	13mm Starter Hole
X1C1N226120	10	Midi	25mm	20mm
X1C1N226125	10	Midi	25mm	25mm
X1C1N326113	10	Midi	25mm	13mm Starter Hole
X1C1N136113	10	Maxi	25mm	13mm Starter Hole
X1C1N136120	10	Maxi	25mm	20mm
X1C1N136125	10	Maxi	25mm	25mm
X1C1N236113	10	Maxi	25mm	13mm Starter Hole
X1C1N236120	10	Maxi	25mm	20mm
X1C1N236125	10	Maxi	25mm	25mm
X1C1N336113	10	Maxi	25mm	13mm Starter Hole
X1C1N126213	10	Midi	35mm	13mm Starter Hole
X1C1N126220	10	Midi	35mm	20mm
X1C1N126225	10	Midi	35mm	25mm
X1C1N126229	10	Midi	35mm	29mm
X1C1N126232	10	Midi	35mm	32mm
X1C1N126235	10	Midi	35mm	35mm
X1C1N226213	10	Midi	35mm	13mm Starter Hole
X1C1N226220	10	Midi	35mm	20mm
X1C1N226225	10	Midi	35mm	25mm
X1C1N226229	10	Midi	35mm	29mm
X1C1N226232	10	Midi	35mm	32mm
X1C1N226235	10	Midi	35mm	35mm
X1C1N326213	10	Midi	35mm	13mm Starter Hole
X1C1N136213	10	Maxi	35mm	13mm Starter Hole
X1C1N136220	10	Maxi	35mm	20mm
X1C1N136225	10	Maxi	35mm	25mm
X1C1N136229	10	Maxi	35mm	29mm
X1C1N136232	10	Maxi	35mm	32mm
X1C1N136235	10	Maxi	35mm	35mm

Product Code	Pack Size	Pouch Size	Plateau Size	Aperture Size
X1C1N236213	10	Maxi	35mm	13mm Starter Hole
X1C1N236220	10	Maxi	35mm	20mm
X1C1N236225	10	Maxi	35mm	25mm
X1C1N236229	10	Maxi	35mm	29mm
X1C1N236232	10	Maxi	35mm	32mm
X1C1N236235	10	Maxi	35mm	35mm
X1C1N336213	10	Maxi	35mm	13mm Starter Hole
X1C1N126313	10	Midi	48mm	13mm Starter Hole
X1C1N126332	10	Midi	48mm	32mm
X1C1N126335	10	Midi	48mm	35mm
X1C1N126338	10	Midi	48mm	38mm
X1C1N126344	10	Midi	48mm	44mm
X1C1N226313	10	Midi	48mm	13mm Starter Hole
X1C1N226332	10	Midi	48mm	32mm
X1C1N226335	10	Midi	48mm	35mm
X1C1N226338	10	Midi	48mm	38mm
X1C1N226344	10	Midi	48mm	44mm
X1C1N326313	10	Midi	48mm	13mm Starter Hole
X1C1N136313	10	Maxi	48mm	13mm Starter Hole
X1C1N136332	10	Maxi	48mm	32mm
X1C1N136335	10	Maxi	48mm	35mm
X1C1N136338	10	Maxi	48mm	38mm
X1C1N136344	10	Maxi	48mm	44mm
X1C1N236313	10	Maxi	48mm	13mm Starter Hole
X1C1N236332	10	Maxi	48mm	32mm
X1C1N236335	10	Maxi	48mm	35mm
X1C1N236338	10	Maxi	48mm	38mm
X1C1N236344	10	Maxi	48mm	44mm
X1C1N336313	10	Maxi	48mm	13mm Starter Hole
X1C1N146413	10	Midi Plus	60mm	13mm Starter Hole
X1C1N246413	10	Midi Plus	60mm	13mm Starter Hole
X1C1N346413	10	Midi Plus	60mm	13mm Starter Hole
X1C1N156413	10	Maxi Plus	60mm	13mm Starter Hole
X1C1N256413	10	Maxi Plus	60mm	13mm Starter Hole
X1C1N356413	10	Maxi Plus	60mm	13mm Starter Hole






Aurum® Plus Soft Convex Drainable

Product Code	Pack Size	Pouch Size	Plateau Size	Aperture Size
X1D1N126113	10	Midi	25mm	13mm Starter Hole
X1D1N126120	10	Midi	25mm	20mm
X1D1N126125	10	Midi	25mm	25mm
X1D1N226113	10	Midi	25mm	13mm Starter Hole
X1D1N226120	10	Midi	25mm	20mm
X1D1N226125	10	Midi	25mm	25mm
X1D1N326113	10	Midi	25mm	13mm Starter Hole
X1D1N136113	10	Maxi	25mm	13mm Starter Hole
X1D1N136120	10	Maxi	25mm	20mm
X1D1N136125	10	Maxi	25mm	25mm
X1D1N236113	10	Maxi	25mm	13mm Starter Hole
X1D1N236120	10	Maxi	25mm	20mm
X1D1N236125	10	Maxi	25mm	25mm
X1D1N336113	10	Maxi	25mm	13mm Starter Hole
X1D1N126213	10	Midi	35mm	13mm Starter Hole
X1D1N126220	10	Midi	35mm	20mm
X1D1N126225	10	Midi	35mm	25mm
X1D1N126229	10	Midi	35mm	29mm
X1D1N126232	10	Midi	35mm	32mm
X1D1N126235	10	Midi	35mm	35mm
X1D1N226213	10	Midi	35mm	13mm Starter Hole
X1D1N226220	10	Midi	35mm	20mm
X1D1N226225	10	Midi	35mm	25mm
X1D1N226229	10	Midi	35mm	29mm
X1D1N226232	10	Midi	35mm	32mm
X1D1N226235	10	Midi	35mm	35mm
X1D1N326213	10	Midi	35mm	13mm Starter Hole
X1D1N136213	10	Maxi	35mm	13mm Starter Hole
X1D1N136220	10	Maxi	35mm	20mm
X1D1N136225	10	Maxi	35mm	25mm
X1D1N136229	10	Maxi	35mm	29mm
X1D1N136232	10	Maxi	35mm	32mm
X1D1N136235	10	Maxi	35mm	35mm

Product Code	Pack Size	Pouch Size	Plateau Size	Aperture Size
X1D1N236213	10	Maxi	35mm	13mm Starter Hole
X1D1N236220	10	Maxi	35mm	20mm
X1D1N236225	10	Maxi	35mm	25mm
X1D1N236229	10	Maxi	35mm	29mm
X1D1N236232	10	Maxi	35mm	32mm
X1D1N236235	10	Maxi	35mm	35mm
X1D1N336213	10	Maxi	35mm	13mm Starter Hole
X1D1N126313	10	Midi	48mm	13mm Starter Hole
X1D1N126332	10	Midi	48mm	32mm
X1D1N126335	10	Midi	48mm	35mm
X1D1N126338	10	Midi	48mm	38mm
X1D1N126344	10	Midi	48mm	44mm
X1D1N226313	10	Midi	48mm	13mm Starter Hole
X1D1N226332	10	Midi	48mm	32mm
X1D1N226335	10	Midi	48mm	35mm
X1D1N226338	10	Midi	48mm	38mm
X1D1N226344	10	Midi	48mm	44mm
X1D1N326313	10	Midi	48mm	13mm Starter Hole
X1D1N136313	10	Maxi	48mm	13mm Starter Hole
X1D1N136332	10	Maxi	48mm	32mm
X1D1N136335	10	Maxi	48mm	35mm
X1D1N136338	10	Maxi	48mm	38mm
X1D1N136344	10	Maxi	48mm	44mm
X1D1N236313	10	Maxi	48mm	13mm Starter Hole
X1D1N236332	10	Maxi	48mm	32mm
X1D1N236335	10	Maxi	48mm	35mm
X1D1N236338	10	Maxi	48mm	38mm
X1D1N236344	10	Maxi	48mm	44mm
X1D1N336313	10	Maxi	48mm	13mm Starter Hole
X1D1N146413	10	Midi Plus	60mm	13mm Starter Hole
X1D1N246413	10	Midi Plus	60mm	13mm Starter Hole
X1D1N346413	10	Midi Plus	60mm	13mm Starter Hole
X1D1N156413	10	Maxi Plus	60mm	13mm Starter Hole
X1D1N256413	10	Maxi Plus	60mm	13mm Starter Hole
X1D1N356413	10	Maxi Plus	60mm	13mm Starter Hole



KEY

-  Black colour pouch
-  Sand colour pouch
-  Clear pouch



Aurum® Plus Soft Convex Urostomy

Product Code	Pack Size	Pouch Size	Plateau Size	Aperture Size
X1U1R126113	10	Midi	25mm	13mm Starter Hole
X1U1R126119	10	Midi	25mm	19mm
X1U1R126122	10	Midi	25mm	22mm
X1U1R126125	10	Midi	25mm	25mm
X1U1R226113	10	Midi	25mm	13mm Starter Hole
X1U1R226119	10	Midi	25mm	19mm
X1U1R226122	10	Midi	25mm	22mm
X1U1R226125	10	Midi	25mm	25mm
X1U1R326113	10	Midi	25mm	13mm Starter Hole
X1U1R136113	10	Maxi	25mm	13mm Starter Hole
X1U1R136119	10	Maxi	25mm	19mm
X1U1R136122	10	Maxi	25mm	22mm
X1U1R136125	10	Maxi	25mm	25mm
X1U1R236113	10	Maxi	25mm	13mm Starter Hole
X1U1R236119	10	Maxi	25mm	19mm
X1U1R236122	10	Maxi	25mm	22mm
X1U1R236125	10	Maxi	25mm	25mm
X1U1R126213	10	Midi	35mm	13mm Starter Hole
X1U1R126219	10	Midi	35mm	19mm
X1U1R126222	10	Midi	35mm	22mm
X1U1R126225	10	Midi	35mm	25mm
X1U1R126229	10	Midi	35mm	29mm
X1U1R126232	10	Midi	35mm	32mm
X1U1R126235	10	Midi	35mm	35mm
X1U1R226213	10	Midi	35mm	13mm Starter Hole
X1U1R226219	10	Midi	35mm	19mm
X1U1R226222	10	Midi	35mm	22mm
X1U1R226225	10	Midi	35mm	25mm
X1U1R226229	10	Midi	35mm	29mm
X1U1R226232	10	Midi	35mm	32mm
X1U1R226235	10	Midi	35mm	35mm
X1U1R326213	10	Midi	35mm	13mm Starter Hole

Product Code	Pack Size	Pouch Size	Plateau Size	Aperture Size
X1U1R136213	10	Maxi	35mm	13mm Starter Hole
X1U1R136219	10	Maxi	35mm	19mm
X1U1R136222	10	Maxi	35mm	22mm
X1U1R136225	10	Maxi	35mm	25mm
X1U1R136229	10	Maxi	35mm	29mm
X1U1R136232	10	Maxi	35mm	32mm
X1U1R136235	10	Maxi	35mm	35mm
X1U1R236213	10	Maxi	35mm	13mm Starter Hole
X1U1R236219	10	Maxi	35mm	19mm
X1U1R236222	10	Maxi	35mm	22mm
X1U1R236225	10	Maxi	35mm	25mm
X1U1R236229	10	Maxi	35mm	29mm
X1U1R236232	10	Maxi	35mm	32mm
X1U1R236235	10	Maxi	35mm	35mm
X1U1R336213	10	Maxi	35mm	13mm Starter Hole
X1U1R126313	10	Midi	48mm	13mm Starter Hole
X1U1R126332	10	Midi	48mm	32mm
X1U1R126335	10	Midi	48mm	35mm
X1U1R226313	10	Midi	48mm	13mm Starter Hole
X1U1R226332	10	Midi	48mm	32mm
X1U1R226335	10	Midi	48mm	35mm
X1U1R326313	10	Midi	48mm	13mm Starter Hole
X1U1R136313	10	Maxi	48mm	13mm Starter Hole
X1U1R136332	10	Maxi	48mm	32mm
X1U1R136335	10	Maxi	48mm	35mm
X1U1R236313	10	Maxi	48mm	13mm Starter Hole
X1U1R236332	10	Maxi	48mm	32mm
X1U1R236335	10	Maxi	48mm	35mm
X1U1R336313	10	Maxi	48mm	13mm Starter Hole






Aurum® Plus Flat Closed

Product Code	Pack Size	Pouch Size	Aperture Size
X1C1A11S	30	Mini	13mm starter hole (cut to 50mm x 70mm)
X1C1A21S	30	Mini	13mm starter hole (cut to 50mm x 70mm)
X1C1A12S	30	Midi	13mm starter hole (cut to 60mm x 80mm)
X1C1A1220	30	Midi	20mm
X1C1A1225	30	Midi	25mm
X1C1A1229	30	Midi	29mm
X1C1A1232	30	Midi	32mm
X1C1A1235	30	Midi	35mm
X1C1A1238	30	Midi	38mm
X1C1A1244	30	Midi	44mm
X1C1A22S	30	Midi	13mm starter hole (cut to 60mm x 80mm)
X1C1A2220	30	Midi	20mm
X1C1A2225	30	Midi	25mm
X1C1A2229	30	Midi	29mm
X1C1A2232	30	Midi	32mm
X1C1A2235	30	Midi	35mm
X1C1A2238	30	Midi	38mm
X1C1A2244	30	Midi	44mm

Product Code	Pack Size	Pouch Size	Aperture Size
X1C1A32S	30	Midi	13mm starter hole (cut to 60mm x 80mm)
X1C1A13S	30	Maxi	13mm starter hole (cut to 60mm x 80mm)
X1C1A1320	30	Maxi	20mm
X1C1A1325	30	Maxi	25mm
X1C1A1329	30	Maxi	29mm
X1C1A1332	30	Maxi	32mm
X1C1A1335	30	Maxi	35mm
X1C1A1338	30	Maxi	38mm
X1C1A1344	30	Maxi	44mm
X1C1A23S	30	Maxi	13mm starter hole (cut to 60mm x 80mm)
X1C1A2320	30	Maxi	20mm
X1C1A2325	30	Maxi	25mm
X1C1A2329	30	Maxi	29mm
X1C1A2332	30	Maxi	32mm
X1C1A2335	30	Maxi	35mm
X1C1A2338	30	Maxi	38mm
X1C1A2344	30	Maxi	44mm
X1C1A33S	30	Maxi	13mm starter hole (cut to 60mm x 80mm)



KEY

-  Black colour pouch
-  Sand colour pouch
-  Clear pouch

Aurum® Plus Flat Drainable

Product Code	Pack Size	Pouch Size	Aperture Size
X1D1A11S	30	Mini	13mm starter hole (cut to 50mm x 70mm)
X1D1A21S	30	Mini	13mm starter hole (cut to 50mm x 70mm)
X1D1A12S	30	Midi	13mm starter hole (cut to 60mm x 80mm)
X1D1A1220	30	Midi	20mm
X1D1A1225	30	Midi	25mm
X1D1A1229	30	Midi	29mm
X1D1A1232	30	Midi	32mm
X1D1A1235	30	Midi	35mm
X1D1A1238	30	Midi	38mm
X1D1A1244	30	Midi	44mm
X1D1A22S	30	Midi	13mm starter hole (cut to 60mm x 80mm)
X1D1A2220	30	Midi	20mm
X1D1A2225	30	Midi	25mm
X1D1A2229	30	Midi	29mm
X1D1A2232	30	Midi	32mm
X1D1A2235	30	Midi	35mm
X1D1A2238	30	Midi	38mm
X1D1A2244	30	Midi	44mm

Product Code	Pack Size	Pouch Size	Aperture Size
X1D1A32S	30	Midi	13mm starter hole (cut to 60mm x 80mm)
X1D1A13S	30	Maxi	13mm starter hole (cut to 60mm x 80mm)
X1D1A1320	30	Maxi	20mm
X1D1A1325	30	Maxi	25mm
X1D1A1329	30	Maxi	29mm
X1D1A1332	30	Maxi	32mm
X1D1A1335	30	Maxi	35mm
X1D1A1338	30	Maxi	38mm
X1D1A1344	30	Maxi	44mm
X1D1A23S	30	Maxi	13mm starter hole (cut to 60mm x 80mm)
X1D1A2320	30	Maxi	20mm
X1D1A2325	30	Maxi	25mm
X1D1A2329	30	Maxi	29mm
X1D1A2332	30	Maxi	32mm
X1D1A2335	30	Maxi	35mm
X1D1A2338	30	Maxi	38mm
X1D1A2344	30	Maxi	44mm
X1D1A33S	30	Maxi	13mm starter hole (cut to 60mm x 80mm)



Aurum® Plus Flat Urostomy

Product Code	Pack Size	Pouch Size	Aperture Size
X1U1E12S	10	Midi	13mm starter hole (cut to 55mm)
X1U1E22S	10	Midi	13mm starter hole (cut to 55mm)
X1U1E32S	10	Midi	13mm starter hole (cut to 55mm)
X1U1E13S	10	Maxi	13mm starter hole (cut to 55mm)
X1U1E1320	10	Maxi	20mm
X1U1E1325	10	Maxi	25mm
X1U1E1329	10	Maxi	29mm
X1U1E1332	10	Maxi	32mm
X1U1E1335	10	Maxi	35mm
X1U1E1338	10	Maxi	38mm
X1U1E23S	10	Maxi	13mm starter hole (cut to 55mm)
X1U1E2320	10	Maxi	20mm
X1U1E2325	10	Maxi	25mm
X1U1E2329	10	Maxi	29mm
X1U1E2332	10	Maxi	32mm
X1U1E2335	10	Maxi	35mm
X1U1E2338	10	Maxi	38mm
X1U1E33S	10	Maxi	13mm starter hole (cut to 55mm)





Order your samples at
www.wellandmedical.com/aurumplus



Welland Medical Ltd., Hydehurst Lane, Crawley, West Sussex RH10 9AS United Kingdom

Tel: +44 (0) 1293 615455 **Fax:** +44 (0) 1293 615411 **Email:** info@wellandmedical.com

©Welland Medical 2023 • Welland Medical Ltd. is a CliniMed® Group Company • Welland®, the Welland logo™ and Aurum® are trademarks of CliniMed (Holdings) Ltd.

www.wellandmedical.com