

Aurum2 
with **Manuka Honey**



Introducing the
Aurum[®]2 range



Innovators in **Stomacare**

Aurum²



with **Manuka Honey**

The new Aurum² range of two-piece colostomy, ileostomy and urostomy pouches and flanges offer new and enhanced features that help improve the quality of lives of ostomates around the world.

One of the key new features of the Aurum² range is the introduction of medical grade Manuka Honey to our skin friendly Hyperflex[®] hydrocolloid, which may help to promote healthy skin around the stoma.

Benefits at a glance...

Confidence

Improved security

The innovative dual security system allows the pouch to be attached both via a flexible rim and adhesive ring

Reduced leakage

Dual-Carb^{®i} (ileostomy) filter, features a new carbon structure and PTFE membrane which repels liquid and oily substances

Comfort

Skin friendly

Welland Medical's tried and tested Hyperflex[®] hydrocolloid with Manuka Honey is kind to even the most vulnerable skin. Natural ingredients help reduce skin irritation

Less ballooning

Improved deodorisation and less ballooning with the new Dual-Carb^{®2} (colostomy) filter

Care

Smooth removal

The easy grip release tab on the pouch facilitates removal of the pouch from the flange

easy2peel[™] tab supports the smooth removal of the flange from the body reducing trauma to the skin

Why Manuka Honey?

Manuka Honey is pure and natural made only from the nectar of the Manuka flower in New Zealand. Historical research has shown this type of active Manuka Honey has characteristics that kill bacteria, suppress inflammation and stimulate the growth of cells, which may aid the healing process.



Colostomy

The new Aurum[®]2 range of closed pouches features the introduction of the split back cover, allowing easy positioning of the pouch and the new Dual-Carb[®]2 filter, specifically designed with colostomists in mind.



Aurum[®]2 Colostomy: Key Benefits

- Skin friendly Hyperflex[®] hydrocolloid with Manuka Honey
- Improved deodorisation and less ballooning with the Dual-Carb[®]2 filter
- A choice of rings and pouch sizes are available with either flat or convex flanges to suit individual needs
- Less skin trauma, due to minimal pressure needed to securely attach the pouch to the flange
- Increased confidence is achieved with the dual security system, whereby the pouch is attached via a flexible connecting rim and adhesive ring
- Easier pouch positioning and output monitoring via the split back cover
- Reduced skin trauma with the easy2peel[™] tab, which aids the smooth removal of the flange
- Quick and simple pouch changing aided by the easy to grip release tab
- Experience exceptional comfort, even during extended wear, with the super soft, quick drying and water repellent cover
- Reduced skin irritation and potential leakage as the cutting guide allows a perfect fit to different stoma shapes and sizes





NEW Dual-Carb[®] 2 filter

How does it work?

The new Dual-Carb[®] 2 filter has an improved carbon formulation which increases odour adsorption for better deodorisation and improved permeability helps to reduce ballooning.



Ileostomy

The new Aurum[®]2 range of drainable pouches features the introduction of the split back cover, allowing easy positioning of the pouch and the new Dual-Carb[®]i filter, specifically designed with ileostomists in mind.



Aurum[®]2 Ileostomy: Key Benefits

- Skin friendly Hyperflex[®] hydrocolloid with Manuka Honey
- Reduced chance of leakage as the Dual-Carb[®]i (ileostomy) filter repels liquid and oily substances
- A choice of rings and pouch sizes are available with either flat or convex flanges to suit individual needs
- Less skin trauma, due to minimal pressure needed to securely attach the pouch to the flange
- Increased confidence is achieved with the dual security system, whereby the pouch is attached via a flexible connecting rim and adhesive ring
- Reduced skin trauma with the easy2peel[™] tab, which aids the smooth removal of the flange
- Easier emptying and cleaning with Welland Medical's innovative integral closure system
- Experience exceptional comfort, even during extended wear, with the super soft, quick drying and water repellent cover
- Easier pouch positioning and output monitoring via the split back cover
- Quick and simple pouch changing aided by the easy to grip release tab
- Reduced skin irritation and potential leakage as the cutting guide allows a perfect fit to different stoma shapes and sizes





NEW Dual-Carb[®] i filter

How does it work?

Dual-Carb[®] i (ileostomy) filter, features a new carbon structure which provides a larger surface area for the carbon, allowing the free flow of air and eliminating odour. A key feature of the filter is the PTFE membrane which repels liquid and oily substances, reducing the chance of leakage.



Urostomy

Have the confidence to wear your pouch and continue your daily activities with the Aurum[®]2 urostomy range.

Offering new and enhanced features, Aurum[®]2 includes Manuka Honey in our skin friendly Hyperflex[®] hydrocolloid and a new convex flange addition – ideal for those with poorly sited, recessed or retracted stomas.



Aurum[®]2 Urostomy: Key Benefits

- Skin friendly Hyperflex[®] hydrocolloid with Manuka Honey
- A choice of rings and pouch sizes are available with either flat or new convex flanges to suit individual needs
- Reduced pouch profile due to the improved tap pocket
- Less skin trauma, due to minimal pressure needed to securely attach the pouch to the flange
- Increased confidence is achieved with the dual security system, whereby the pouch is attached via a flexible connecting rim and adhesive ring
- Reduced risk of infection as the non-return valve prevents the reflux of urine
- Easier pouch positioning and output monitoring via the split back cover
- Reduced skin trauma with the easy2peel[™] tab, which aids the smooth removal of the flange
- Quick and simple pouch changing aided by the easy to grip release tab
- Experience exceptional comfort, even during extended wear, with the super soft, quick drying and water repellent cover
- Increased comfort with the extremely soft and flexible tap



NEW
Convex
flange for
urostomy



Aurum2 
with **Manuka Honey**

Product Codes

Aurum[®]2 Two Piece Colostomy Pouches

| Product Code | Pack Size | Pouch Size | Ring Size | Beige or Clear |
|--------------|-----------|------------|-----------|----------------|
| XMH2C445 | 30 | Mini | 45mm | Beige |
| XMH2C545 | 30 | Midi | 45mm | Beige |
| XMH2C145 | 30 | Maxi | 45mm | Beige |
| XMH2C455 | 30 | Mini | 55mm | Beige |
| XMH2C555 | 30 | Midi | 55mm | Beige |
| XMH2C155 | 30 | Maxi | 55mm | Beige |
| XMH2C570 | 10 | Midi | 70mm | Beige |
| XMH2C170 | 10 | Maxi | 70mm | Beige |



Aurum[®]2 Two Piece Ileostomy Pouches

| Product Code | Pack Size | Pouch Size | Ring Size | Beige or Clear |
|--------------|-----------|------------|-----------|----------------|
| XMH2D445 | 30 | Mini | 45mm | Beige |
| XMH2D545 | 30 | Midi | 45mm | Beige |
| XMH2D145 | 30 | Maxi | 45mm | Beige |
| XMH2D455 | 30 | Mini | 55mm | Beige |
| XMH2D555 | 30 | Midi | 55mm | Beige |
| XMH2D155 | 30 | Maxi | 55mm | Beige |
| XMH2D570 | 10 | Midi | 70mm | Beige |
| XMH2D170 | 10 | Maxi | 70mm | Beige |



Aurum[®]2 Two Piece Urostomy Pouches

| Product Code | Pack Size | Pouch Size | Ring Size | Beige or Clear |
|--------------|-----------|------------|-----------|----------------|
| XMH2U545 | 10 | Midi | 45mm | Beige |
| XMH2U145 | 10 | Maxi | 45mm | Beige |
| XMH2U555 | 10 | Midi | 55mm | Beige |
| XMH2U155 | 10 | Maxi | 55mm | Beige |



Aurum² Two Piece Flat Flanges

| Product Code | Pack Size | Ring Size | Aperture Size |
|--------------|-----------|-----------|---------------|
| XMH2F413 | 5 | 45mm | 13mm - 40mm |
| XMH2F419 | 5 | 45mm | 19mm |
| XMH2F422 | 5 | 45mm | 22mm |
| XMH2F425 | 5 | 45mm | 25mm |
| XMH2F429 | 5 | 45mm | 29mm |
| XMH2F432 | 5 | 45mm | 32mm |
| XMH2F435 | 5 | 45mm | 35mm |
| XMH2F438 | 5 | 45mm | 38mm |
| XMH2F513 | 5 | 55mm | 13mm - 50mm |
| XMH2F538 | 5 | 55mm | 38mm |
| XMH2F544 | 5 | 55mm | 44mm |
| XMH2F713 | 5 | 70mm | 13mm - 65mm |



Aurum² Two Piece Convex Flanges

| Product Code | Pack Size | Ring Size | Aperture Size |
|--------------|-----------|-----------|---------------|
| XMH2NF413 | 5 | 45mm | 13mm - 25mm |
| XMH2NF419 | 5 | 45mm | 19mm |
| XMH2NF422 | 5 | 45mm | 22mm |
| XMH2NF513 | 5 | 55mm | 13mm - 35mm |
| XMH2NF525 | 5 | 55mm | 25mm |
| XMH2NF529 | 5 | 55mm | 29mm |
| XMH2NF532 | 5 | 55mm | 32mm |
| XMH2NF535 | 5 | 55mm | 35mm |
| XMH2NF713 | 5 | 70mm | 13mm - 48mm |



Aurum² Two Piece Urostomy Flat Flanges

| Product Code | Pack Size | Ring Size | Aperture Size |
|--------------|-----------|-----------|---------------|
| XMH2UF413 | 5 | 45mm | 13mm - 40mm |
| XMH2UF419 | 5 | 45mm | 19mm |
| XMH2UF422 | 5 | 45mm | 22mm |
| XMH2UF425 | 5 | 45mm | 25mm |
| XMH2UF429 | 5 | 45mm | 29mm |
| XMH2UF432 | 5 | 45mm | 32mm |
| XMH2UF435 | 5 | 45mm | 35mm |
| XMH2UF438 | 5 | 45mm | 38mm |
| XMH2UF513 | 5 | 55mm | 13mm - 50mm |

Aurum² Two Piece Urostomy Convex Flanges

| Product Code | Pack Size | Ring Size | Aperture Size |
|--------------|-----------|-----------|---------------|
| XMH2NUF413 | 5 | 45mm | 13mm - 25mm |
| XMH2NUF419 | 5 | 45mm | 19mm |
| XMH2NUF422 | 5 | 45mm | 22mm |
| XMH2NUF513 | 5 | 55mm | 13mm - 35mm |
| XMH2NUF525 | 5 | 55mm | 25mm |
| XMH2NUF529 | 5 | 55mm | 29mm |
| XMH2NUF532 | 5 | 55mm | 32mm |
| XMH2NUF535 | 5 | 55mm | 35mm |

Aurum2

with **Manuka Honey**



Welland Medical Ltd., Hydehurst Lane, Crawley, West Sussex RH10 9AS United Kingdom
Tel: +44 (0) 1293 615455 **Fax:** +44 (0) 1293 615411 **Email:** info@wellandmedical.com (A CliniMed® Group Company)
All products marked ®, Welland® and the Welland logo™ are trade marks of CliniMed (Holdings) Limited.

www.wellandmedical.com

