

Profiled to your body

- The surface is smooth without any hard, rigid edges that may potentially harm the stoma
- Can be folded during application ensuring a secure fit across the contours of the body
- Conforms gently around or into skin folds, scar tissue, moats or recesses at body temperature, reducing the need for pastes or washers
- The Hyperflex® hydrocolloid has added Manuka honey, which may help to promote healthy skin around the stoma
- Can be inverted to shape over a parastomal hernia



Aurum Profile
with Manuka Honey

Aurum Profile
with Manuka Honey

Introducing the Aurum® Profile range



Aurum Profile
with Manuka Honey



Welland Medical Ltd., Hydehurst Lane, Crawley, West Sussex RH10 9AS United Kingdom
Tel: +44 (0) 1293 615455 Fax: +44 (0) 1293 615411 Email: info@wellandmedical.com (A CliniMed® Group Company)
All products marked "Welland" and the Welland logo™ are trade marks of CliniMed (Holdings) Limited.

www.wellandmedical.com



PID 4696



Innovators in Stomacare



Soft and mouldable,
adapts to any body profile.

Ostomates with uneven peristomal skin can often experience leaks and sore skin around the stoma when using a flat flange.

Aurum® Profile has a soft and mouldable plateau that provides a secure seal around these complex body contours where a flat flange is not suitable. This reduces the need for paste or washers; giving the ostomate confidence to have a secure, leak-free fit.

Benefits at a glance...

Skin Friendly

The Hyperflex® hydrocolloid flange has added Manuka honey, which may help to promote healthy skin around the stoma.

Flexible Solution

The soft and mouldable flange reduces the risk of discomfort and skin trauma.

Manuka Honey

This type of active Manuka honey has characteristics that kill bacteria, suppress inflammation and stimulate the growth of cells.

Ideal For

- Parastomal hernias
- Peristomal irregularities

Ileostomy Product Codes



Aurum® Profile Drainable Pouches

Product Code	Pack Size	Plateau Size	Pouch Size	Aperture Size	Beige or Clear
XMHCDS413	10	Small	Mini	13mm – 25mm	Beige
XMHCDS613	10	Small	Mini	13mm – 25mm	Clear
XMHCDS513	10	Small	Midi	13mm – 25mm	Beige
XMHCDS519	10	Small	Midi	19mm	Beige
XMHCDS522	10	Small	Midi	22mm	Beige
XMHCDS525	10	Small	Midi	25mm	Beige
XMHCDS713	10	Small	Midi	13mm – 25mm	Clear
XMHCDS913	10	Small	Maxi	13mm – 25mm	Beige
XMHCDS919	10	Small	Maxi	19mm	Beige
XMHCDS922	10	Small	Maxi	22mm	Beige
XMHCDS925	10	Small	Maxi	25mm	Beige
XMHCDS313	10	Small	Maxi	13mm – 25mm	Clear
XMHCDS413	10	Medium	Mini	13mm – 35mm	Beige
XMHCDS613	10	Medium	Mini	13mm – 35mm	Clear
XMHCDS513	10	Medium	Midi	13mm – 35mm	Beige
XMHCDS519	10	Medium	Midi	19mm	Beige
XMHCDS522	10	Medium	Midi	22mm	Beige
XMHCDS525	10	Medium	Midi	25mm	Beige
XMHCDS529	10	Medium	Midi	29mm	Beige
XMHCDS532	10	Medium	Midi	32mm	Beige
XMHCDS535	10	Medium	Midi	35mm	Beige
XMHCDS713	10	Medium	Midi	13mm – 35mm	Clear
XMHCDS913	10	Medium	Maxi	13mm – 35mm	Beige
XMHCDS919	10	Medium	Maxi	19mm	Beige

Urostomy Product Codes



Aurum® Profile Urostomy Pouches

Product Code	Pack Size	Plateau Size	Pouch Size	Aperture Size	Beige or Clear
XMHCUS513	10	Small	Midi	13mm – 25mm	Beige
XMHCUS519	10	Small	Midi	19mm	Beige
XMHCUS522	10	Small	Midi	22mm	Beige
XMHCUS525	10	Small	Midi	25mm	Beige
XMHCUS713	10	Small	Midi	13mm – 25mm	Clear
XMHCUS913	10	Small	Maxi	13mm – 25mm	Beige
XMHCUS919	10	Small	Maxi	19mm	Beige
XMHCUS922	10	Small	Maxi	22mm	Beige
XMHCUS925	10	Small	Maxi	25mm	Beige
XMHCUS313	10	Small	Maxi	13mm – 25mm	Clear
XMHCUS513	10	Medium	Midi	13mm – 35mm	Beige
XMHCUS519	10	Medium	Midi	19mm	Beige
XMHCUS522	10	Medium	Midi	22mm	Beige
XMHCUS525	10	Medium	Midi	25mm	Beige
XMHCUS529	10	Medium	Midi	29mm	Beige
XMHCUS532	10	Medium	Midi	32mm	Beige
XMHCUS535	10	Medium	Midi	35mm	Beige
XMHCUM713	10	Medium	Midi	13mm – 35mm	Clear
XMHCUM913	10	Medium	Maxi	13mm – 35mm	Beige
XMHCUM919	10	Medium	Maxi	19mm	Beige

Product Code	Pack Size	Plateau Size	Pouch Size	Aperture Size	Beige or Clear
XMHCUM922	10	Medium	Maxi	22mm	Beige
XMHCUM925	10	Medium	Maxi	25mm	Beige
XMHCUM929	10	Medium	Maxi	29mm	Beige
XMHCUM932	10	Medium	Maxi	32mm	Beige
XMHCUM935	10	Medium	Maxi	35mm	Beige
XMHCUM313	10	Medium	Maxi	13mm – 35mm	Clear
XMHCUL513	10	Large	Midi	13mm – 48mm	Beige
XMHCUL532	10	Large	Midi	32mm	Beige
XMHCUL535	10	Large	Midi	35mm	Beige
XMHCUL538	10	Large	Midi	38mm	Beige
XMHCUL544	10	Large	Midi	44mm	Beige
XMHCUL713	10	Large	Midi	13mm – 48mm	Clear
XMHCUL913	10	Large	Maxi	13mm – 48mm	Beige
XMHCUL932	10	Large	Maxi	32mm	Beige
XMHCUL935	10	Large	Maxi	35mm	Beige
XMHCUL938	10	Large	Maxi	38mm	Beige
XMHCUL313	10	Large	Maxi	13mm – 48mm	Clear

The Flexible Flange